



Sterile Eye Wash Solution 500ml

SKU: REL-900

[\(View Product Online\)](#)

Description:

The Eye Wash Solution 500ml is designed for the proper cleansing and irrigation of the eyes. Each bottle contains 500ml of sterile solution, which is ideal for treating minor eye injuries.

Mandatory under HSE regulations.

Sterile, safe, and non-toxic.

Twist seal cap.

- Sterile, 500ml bottle of eyewash solution, Sodium chloride 0.9%

Frequently Asked Questions:

Q: What is an eye wash solution for?

A: This sterile eye wash is made for cleansing and irrigation of the eyes; removing any contaminants like dust or debris. It aids in treating minor eye injuries.

Q: How should the eye wash solution be used in an emergency?

A: In the event of minor eye injury, or eye contamination, position the person with their slightly head tilted back. Open the sterile eye wash, and aim the nozzle at the inner corner of the affected eye, and gently squeeze the bottle of eye wash so that the solution flows across the eye, flushing out contaminants. Continue for at least 15 minutes or until professional medical help is available. Ensure to gently hold the eyelids of the affected eye open, and ask the person to keep moving their eyeball to ensure thorough rinsing with the eye wash solution.

Q: What is the shelf life of the eye wash solution?

A: The eye wash solution typically has a shelf life of up to 4-5 years.

Q: Are there any precautions to consider when using this eye wash solution?

A: Make sure the bottle of eye wash solution stays sealed until it is ready to be used to keep the best sterility. Do not use the eye wash solution if the seal is broken or the solution is discoloured.

Q: Can this sterile eye wash solution be used for other purposes?

A: This sterile eye wash solution is specifically made for eye irrigation and minor eye injuries. For wound cleansing or other first aid needs, consider our other first aid supplies such as our cotton wool for gentle cleaning.

Q: How should the eye wash solution be stored?

A: It is recommended to store the bottle of sterile eye wash in a cool, dry place away from direct sunlight or any extreme temperatures to preserve sterility and optimal effectiveness.

Q: Is any training required to use eye wash solution?

A: While the eye wash solution is very user-friendly, basic first aid training can improve the familiarity of the product as well the effectiveness of its use. We also recommend the user-friendly eye wash kit, stocked with everything you would need for cleansing and irrigation of the eyes.